


Non-invasive and high-throughput electrophysiological analysis of apical and basolateral membranes of RPE cells

Colby F. Lewallen, Ph.D.
National Eye Institute
PI: Dr. Kapil Bharti

Outline


- Background of RPE electrophysiology
- Derivation of impedance model constraint
- Validation of impedance model constraint
- Testing using LCA16 iPSC-RPE line

Electrophysiology background



Key Features of RPE

Functions




Characteristics

1. Polarized
2. Pigmented
3. Hexagonal shape
4. Tight junctions
5. Tightly adhered
6. Non-proliferative



Validation of iPSC-Derived RPE



Characteristics Validation

- **Cell purity:** flow cytometry
- **Morphology:** quantitative shape metrics
- **Molecular markers:** RPE-specific gene expression
- **Membrane protein localization:** Polarized distribution of proteins (Na⁺/K⁺ ATPase, ezrin, ZO-1)

Functional Validation

- **Polarized secretion:** VEGF secretion assays
- **Phagocytosis & lysosomal function:** Uptake and degradation of photoreceptor outer segments
- **Barrier integrity:** Trans-epithelial resistance
- **Ion transport:** Calcium signaling, ion flux assays (e.g., patch-clamp)
- **Metabolism:** Mitochondrial activity profiling

**Snapshot of
cell quality
and function**

Sharma et al. (2019)
Science Translational Medicine



Validation of iPSC-Derived RPE

Characteristics Validation


- **Cell purity:** flow cytometry
- **Morphology:** quantitative shape metrics
Molecular markers: RPE-specific gene expression
- **Membrane protein localization:** Polarized distribution of proteins (Na⁺/K⁺ ATPase, ezrin, ZO-1)

Functional Validation

- **Polarized secretion:** VEGF secretion assays
- **Phagocytosis & lysosomal function:** Uptake and degradation of photoreceptor outer segments
- **Barrier integrity:** Trans-epithelial resistance
- **Ion transport:** Calcium signaling, ion flux assays (e.g., patch-clamp)
- **Metabolism:** Mitochondrial activity profiling




*High throughput
assessment of these
properties is **expensive**
and **impractical** for
most studies*



How can we **minimize time** and **maximize measured quality indicators** for high throughput analysis?

How can we increase the **sensitivity** of these measurements?

RPE polar properties represented by equivalent circuit




Example image of setup for measurements (historical)


$R_{a/b/s/sol}$ - apical/basal/shunt/solution **resistance**

$C_{a/b}$ - apical/ basal **capacitance**

ATP response of RPE




ATP response of RPE



One assay reveals complex functional health of epithelia



- Epithelial polarization
- Tight junction formation

fc-iPSC-RPE1-3





- Cell-to-cell homogeneity
- Intracellular ionic composition
- General molecular mechanisms



fC-iPSC-RPE1-3





fC-iPSC-RPE1-3



fC-iPSC-RPE1-3





fC-iPSC-RPE1-3




Limitations

1. Difficult (requires training and precision)
2. Low throughput (1-3 measurements/day)
3. Expensive (>\$100k)
4. Unreliable (can lose cells, hard to make pipettes, etc.)
5. Invasive (requires puncturing/touching the cells)

Epithelial polarity facilitates essential retinal function



Current method to measure parameters of epithelial polarity



INVASIVE and LOW-THROUGHPUT


Goal: Develop a **non-invasive** and **high-throughput** method for the electrophysiological analysis of polarized epithelial cells using extracellular electrochemical impedance spectroscopy (EIS)

Derivation of impedance model

Derivation of equivalent impedance $Z_{eq}(w)$


$$Z_a(\omega) = \frac{R_a}{1 + i\omega R_a C_a}$$

$$Z_b(\omega) = \frac{R_b}{1 + i\omega R_b C_b}$$



Imaginary number: $i = \sqrt{-1}$

Derivation of equivalent impedance model: $Z_{eq}(\omega)$



Combining the
circuit elements



$$Z_{eq}(\omega) = R_{blank} + \frac{R_s(Z_a(\omega) + Z_b(\omega))}{R_s + Z_a(\omega) + Z_b(\omega)}$$

Epithelial resistance and capacity parameters can be extrapolated from measured impedance



$$Z_{eq}(\omega) = R_{blank} + \frac{R_s(Z_a(\omega) + Z_b(\omega))}{R_s + Z_a(\omega) + Z_b(\omega)}$$



$$R_{blank} = Z_{eq}(\infty)$$


$$TER = Z_{eq}(0) = R_{blank} + \frac{R_s(R_a + R_b)}{R_s + R_a + R_b}$$

$$\frac{1}{TEC} = \frac{2}{\pi} \int_0^\infty \Re(Z_{eq}(\omega) - R_{blank}) d\omega$$





Resolving underdetermined fit

Principle of fitting




★ Location of best fit

No clear minimum when fitting $Z_{eq}(w)$




Set target C_b during fitting to get distinct solution


- C_b should not change much between cell lines.
- Primarily growth substrate dependent.



Arvydas Maminishkis,
MD, PhD




Simple constraint converges to distinct solution




Validation

Validate constraint accuracy with **SIMULATED** data



simulate



fit



validate

$$Z_{eq}(\omega) = R_{blank} + \frac{R_s(Z_a(\omega) + Z_b(\omega))}{R_s + Z_a(\omega) + Z_b(\omega)}$$

& $C_b = \text{constant}$

Keep the circuit diagram in all pertinent slides

Accurate prediction of membrane parameters is not impacted by accuracy of constrained basal capacitance value



To do: Put sentence explaining the impact of C_b error on the prediction of other membrane parameters

Biology validation: polarized mutation

Test model with biological data (RPE)

KCNJ13 variants cause LCA16, disrupting Kir7.1 and apical RPE K⁺ conductance

Insert schematic of iPSC->RPE



David Gamm, MD, PhD



Bikash Pattnaik, PhD



Omar Memon



Jair Montford



Samuel Ramirez



Casey Cargill

Kir maturation into RPE

Schematic depicting weekly measurements




Casey Cargill


SAM.png

Samuel Ramirez


Asymptomatic family member (heterozygote)



LCA16 Patient (homozygote)



Diminished Ra response to Ba++ indicates lack of functional Kir7.1 in LCA16 patient line.



Casey Cargill³⁵


Conclusion: We generated a **non-invasive** alternative to intracellular electrophysiological assays for **high-throughput** analysis of polarized epithelial cells

Acknowledgements!

To do: add
lab/collaborator
names



Comparison of RtArea across Cell Lines



Questions